**Schlauch-pumpe**

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**REGLO Digital****2 Kanäle**

ISM831 6 Rollen

ISM832 8 Rollen

ISM596 12 Rollen

**4 Kanäle**

ISM833 6 Rollen

ISM834 8 Rollen

ISM597 12 Rollen

mit Settings-Menü

**Tubing pump**

Microprocessor  
controlled

**REGLO Digital****2 channels**

ISM831 6 rollers

ISM832 8 rollers

ISM596 12 rollers

**4 channels**

ISM833 6 rollers

ISM834 8 rollers

ISM597 12 rollers

with menu settings

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**REGLO Digital****2 canaux**

ISM831 6 galets

ISM832 8 galets

ISM596 12 galets


**4 canaux**

ISM833 6 galets

ISM834 8 galets

ISM597 12 galets

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Inhaltsverzeichnis		Contents		Sommaire	
Sicherheitsvorkehrungen	4	Safety precautions	4	Mesures de sécurité	4
Garantiebestimmungen	6	Warranty terms	6	Conditions de garantie	6
Produkt	7	Product	7	Produit	7
Geräterückwand	8	Real panel	8	Panneau arrière	8
Netzspannung	8	Mains voltage	8	Tension d'alimentation	8
Sicherungen wechseln	8	Changing the fuses	8	Remplacement des fusibles	8
Bedienungspanel	9	Operating panel	9	Tableau de commande	9
Start-Information	10	Start-up information	10	Information de mise en route	10
Grundeinstellungen		Basic settings		Réglage de base	
Steuertasten	11	Control keys	11	Touches de commande	11
Wählen	12	Selecting	12	Choisir	12
Cycles (Dosierzyklen)	13	Cycles (Dispensing cycles)	13	Cycles (Cycles de dosage)	13
Backsteps	13	Backsteps	13	Backsteps	13
Adresse (RS232)	13	Address (RS232)	13	Adresse (RS232)	13
Fußschalter	13	Footswitch	13	Pédale de commande	13
Autostart	14	Autostart	14	Autostart	14
Digits (Anzeigen)	14	Digits (Display)	14	Digits (Affichage)	14
Display Speed	14	Display Speed	14	Display Speed	14
Schlauch Innendurchmesser	15	Tube inner diameter	15	Diamètre intérieure du tube	15
Pumpenkopf	15	Pump-head	15	Tête de pompe	15
Initialisierung	16	Initializing	16	Initialisation	16
Schläuche einlegen	17	Inserting the tubing	17	Insertion des tubes	17
Pumpen nach Fließrate	18	Pumping by flow rate	18	Pompage selon le débit	18
Total-volumen	18	Total volume	18	Volume total	18
Fließrate kalibrieren	19	Calibrating the flow rate	19	Calibration du débit	19
Dosieren nach Volumen	20	Dispensing by volume	20	Dosage selon le volume	20
Volumen kalibrieren	21	Calibrating the volume	21	Calibration du volume	21
Default-Kalibration		Default calibration of		Calibration par défaut	
Fließrate	22	flow rate	22	du débit	22
Volumen	22	volume	22	du volume	22
Dosieren nach Zeit	23	Dispensing by time	23	Dosage selon le temps	23

## Hinweis (Kassetten)

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## Please note (cassettes)

When using new tubing for the first time, it may occur that, depending on the tubing used (hardness and diameter), the pump cannot be primed and, hence, does not deliver the liquid. If that is the case, it is recommended that the tubing be wet and the pump run with the tubing inserted for about 15 to 30 minutes. For applications with a high differential pressure (>1 bar), we recommend the use of cassettes with pressure lever.

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## Please note

We recommend you read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

**ISMATEC® does not take liability for any damage resulting from the use of an ISMATEC® pump.**

ISMATEC® does not admit responsibility for the handling of chemicals.

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## Safety precautions

ISMATEC® tubing pumps are designed for pumping and dispensing applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed.

- ▶ The pump must not be operated outside the designed operating and environmental conditions. In particular, it is the user's responsibility not to exceed the specified maximum differential pressure for the pump-head.
- ▶ The pump must not be used:
  - for medical applications on human beings.
  - in explosion proof chambers or in the presence of flammable gases or fumes.
- ▶ The pump must be switched OFF when pump-heads, cassettes or tubing are inserted or changed.
- ▶ The permeability of tubing depends on the material used and pressure conditions. Tubing can also become electro-statically charged. Please be aware of possible hazards when routing tubing in explosion-proof chambers.
- ▶ Pump-heads consist of rotating parts. Therefore, the pump must not be operated before the cassettes are fully snapped-in.

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## Caution

Rotating pumphead creates a pinch and crush hazard.

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## Safety precautions

- ▶ Do not manipulate the pump-head before the pump is switched off and disconnected from the mains supply.
- ▶ Be particularly cautious that no parts of your body such as fingers, long hair, etc. or jewellery, or loose objects such as cables or tubing, etc. can be trapped by the rotating roller head.
- ▶ Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken. When pumping corrosive media we recommend to place the pump in a collecting basin.

⚠ The pump should not come in contact with the liquid.

- ▶ Only new fuses, according to the specifications stated on Page 8 in this manual, must be used.
- ▶ The fuse holder must not be short-circuited.
- ▶ Do not open or remove the housing while the pump is operating.
- ▶ Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- ▶ For service and repairs carried out by the customer or by third-party companies ISMATEC® denies any responsibility.

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### Warranty

For all parts manufactured by ISMATEC®	2 years
All other parts, excluding consumables:	1 year
In case of any queries, please contact your ISMATEC® representative.	

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
Warranty terms
We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.
If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.
A defective pump must be returned in the original ISMATEC® packaging or in a packaging of equal quality.
The duration of the warranty is not affected by making a claim for warranty service.
Further claims are excluded. Shipping costs are charged to the customer.
<b>Our warranty becomes invalid in the case of:</b> <ul style="list-style-type: none"> <li>▶ improper operation by the user, or if the pump is diverted from its proper use.</li> <li>▶ unauthorized modification or misuse by the user or by a third-party.</li> <li>▶ improper site preparation and maintenance.</li> <li>▶ operation outside of the environmental and electrical specifications for the product.</li> <li>▶ use of third-party software, hardware, accessories or consumables purchased by the user and which do not comply with our specifications.</li> </ul>

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#### Contents of the package

- ▶ Pump REGLO *Digital* (type as ordered)
- ▶ 2 or 4 cassettes
- ▶ 1 Power Cord with country specific mains plug
- ▶ Operating Manual

**Please check the package** and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

Complaints can only be accepted within 8 days from receipt of the goods.

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

## Rear panel

1. RS232 IN (female) for operation controlling via PC
  2. RS232 OUT (male) for cascading up to 8 pumps
- Digital input (TTL level)  
Run/Stop
- Digital output 0/10 kHz  
(speed output), Busy (TTL Level)
3. Mains socket
  4. Fuse holder

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## Mains voltage

Mains voltage	Fuse rating
100–240 V <sub>AC</sub> 50/60 Hz	1 x T2A/250V

### ⚠ Socket/Power cord

Use exclusively the originally applied power cord. The socket must be connected to earth ground (protective conductor contact).

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## Changing the fuses

➡ Switch the pump OFF, pull out the mains plug.


1. Pull out the fuse holder by opening it with a small screw-driver (size 0).
2. Remove fuse (a) and replace it with the spare fuse (b) 1 x T2A/250V.

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## Operating panel

(\* = Settings Page 11)

1. On/Off switch
2. Digital LED display
3. Increase value
4. Reduce value  
\* settings
5. RUN/STOP (double function)
  - a. Starts and stops the pump
  - b. Interrupts and continues a dispensing cycle (Page 20)

➡ The function Start/Stop can also be triggered by an external device such as a footswitch. (Accessories on Page 40).

### 6. RESET (double function)

#### a. Rotation direction

Changes the rotation direction. By changing the rotation direction or starting the pump, the rotation direction is indicated in the display with a running hyphen.

#### b. RESET

- Resets the accumulated volume in mode TOTAL to zero (Page 18)
- Interrupts dispensing cycle for good (Page 20)  
\* escape

### 7. CAL/MAX (double function)

- a. CAL** = Calibrating key for the flow rate or dispensing volume
- b. MAX** = Max. speed in the modes PUMP•Flow rate and TOTAL (keep the MAX key pressed down while the pump is running)  
\* ok

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
Selecting operating mode  
\* default

### 9. LEDs for active operating mode

- ▶ **PUMP • Flow rate**  
Pumping by flow rate (ml/min)
- ▶ **DISP • Volume**  
Dispensing by volume (ml)
- ▶ **DISP • Time**  
Dispensing by time
- ▶ **PAUSE • Time**  
Pause time for dispensing  
at intervals
- ▶ **TOTAL**  
Pumping with Read-out of totally  
delivered volume

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## Start-up information

After switching on the power supply  
switch, the following values  
are displayed:

1. LED test »8.8.8.8.«
2. Firmware version
3. Displays the operating mode last used

⚠ **Before starting the pump for  
the first time or after connecting a  
new pump-head**, the ID-code of the  
pump-head must be entered (see basic  
settings, Page 15).

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## Please note

The settings menu requires a break of **5 seconds** between power off and on.

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## Selecting the basic settings

1. Switch the pump OFF and wait for 5 seconds.
  2. Keep the settings key pressed and switch the pump ON.
  3. Select the basic settings by using the ▲▼ keys
    - Cycles (dispensing cycles)
    - Back-steps (number of back-steps for volume dispensing)
    - Address (RS232, pump address)
    - Footswitch (Toggle/Direct)
    - Autostart (Standby/RUN)
    - Digits (Number of digits displayed for flow rate and dispensing volume)
    - Display Speed (Type of speed display)
    - Tube (Inner diameter)
    - Head = ID-Code (pump-head)
    - Init (resets basic settings to the default values except pump-head and tube diameter)
  4. ok = confirm your selection
  5. Select the value by using the ▲▼ keys, or
    - press escape for cancelling wrong input
    - the values can be reset to the default value by using the default key
  6. ok = confirm entered value
  7. escape  
Quit menu basic settings
- ➡ **Basic settings refer to the following pages**

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Number of dispensing steps 0..9999  
(0 = Default value, infinity)

### Backsteps

Number of roller back-steps (for drip-free dispensing in the mode DISP Volume) from 0 .. 100 (0 = Default value)

### Address (RS232)

Pump address 1 – 8 (Page 31)  
(1=default value)

### Footswitch

- ▶ Toggle (default value)
  - Step on footswitch  
= pump starts running
  - Press footswitch again  
= pump stops
- ▶ Direct  
the pump runs as long as the footswitch is pressed

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## Basic settings

### Configuring Autostart

Autostart determines the pump action when the mains supply is on again, i.e. after power black-out or when the pump is controlled via mains power.

#### Standby

Pump is stand-by when switched on (default value)

#### RUN

Pump runs when switched on

### Digits (display)

Number of significant digits, 3 or 4 (3 = default value) when displaying flow rate and dispensing volume.

### Display Speed

Display of the pump speed in mode Pump•Flow rate as

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### Tubing inner diameter

Configuring the tubing inner diameter  
(0.13 mm = default)

**Head** Identification code (ID) of the pump-head, see also table on Page 39 (2-06 = default value)

➡ The pump-head ID is composed of:

- ▶ number of channels (2 or 4)
- ▶ number of rollers (6, 8 or 12)

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## Basic settings

### Init (Initializing)

Pressing the ok key resets the following basic settings to the default values (except pump-head and tube diameter):

Mode	PUMP•Flow rate
Speed	40 rpm
Dispensing volume	30 roller steps
Dispensing time	4.50 sec
Pause time	2.00 sec
Rotation direction	clockwise
Back-steps	0
Number of dispensing cycles	0
Pump-head digits display	3
Footswitch	Toggle
Autostart	Standby
Display	Flow rate

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

### 3-stop color-coded tubing

Please refer to our detailed tubing documentation.

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## Inserting the tubing

⚠ First enter the tube I.D. as described on Page 15.

- ➡ Switch the pump OFF
- Remove the cassette by slightly pressing the fixing-tongue and lifting it simultaneously.
- Insert the 3-stop color-coded tubing with one stopper into the cassette.
- Let the tubing hang down (prevent it from being twisted).
- Insert the tubing with the second stopper at the other end of the cassette.
- Reinsert the cassette into the roller-head.
- ➡ Switch the pump ON

➡ When the pump is idle, release all cassettes at the fixing-tongue (2)

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## Calibrating the flow rate

1. Change mode to PUMP • Flow rate

Enter the required flow rate with the ▲▼ keys

➡ Check whether the correct tubing I.D. is selected (see Page 15)

2. Change mode to DISP • Time

- Enter 60 seconds by using the ▲▼ keys
- Start pump with RUN/STOP

➡ Pump stops automatically after 60 seconds

➡ Measure the dispensed liquid by volume or weight

Calibrate the ascertained value as follows:

3. Return to mode PUMP • Flow rate

4. Press the CAL/MAX-button (displayed value blinks)  
Enter the measured value (balance) with the ▲▼ keys

5. Confirm with the CAL/MAX-key (the flow rate setting is readjusted to the initially entered set point as close as possible). The rotation speed has a resolution of 1/100 rpm. After calibration, the initially set flow rate can only be approached.

➡ Setting the default calibration: see Page 22

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➡ For dispensing accurately, it is recommended that the pump be calibrated (see Page 21).

3. Start pump with RUN/STOP. The pump reduces the rotation speed shortly before the end of the dispensing cycle providing controllable and drop-precise dispensing cycles.

➡ The dispensing speed can be adjusted in the mode PUMP • Flow rate.

➡ With the ▲▼ keys the dispensing volume can also be changed **even during the dispensing process**. The newly entered set point appears briefly in the display.

➡ During the dispensing process the flow rate (PUMP • Flow rate) and the total volume can be retrieved by pressing the MODE key.

➡ The flow rate can be changed even during the dispensing process by using the ▲▼ keys.

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## Calibrating the volume

1. Change the MODE to DISP • Volume  
➡ Use the ▲▼ keys for entering the required dispensing volume
2. Start with RUN/STOP  
➡ The pump reduces the rotation speed shortly before the end of the dispensing cycle providing controllable and drop-precise dispensing volumes.  
➡ Measure the dispensed liquid by volume or weight.
3. Stay in the same mode (DISP • Volume) and press the CAL/MAX-key (display blinks)
4. Enter the effectively dispensed volume with the ▲▼ keys
5. Save with the CAL/MAX key  
➡ The calibrated volume is adjusted to the initially entered set point as close as possible. (The volume per roller step determines the smallest calibration unit. A dispensing volume is always the multiple of this unit. The initially set volume can often only be approached).
6. Start pump with RUN/STOP  
➡ Depending on the application conditions repeated calibration may be necessary.  
➡ The smallest possible dispensing volume corresponds to the volume of one roller step.

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24h

Please note

A dispensing cycle can be interrupted by pressing the RUN/STOP key.

– the volume dispensed up to now is blinking

– by pressing the RUN/STOP key again, the dispensing cycle is continued

– by pressing the RESET-key the dispensing cycle is stopped for good

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Dispensing by time

The dispensing time can be entered from 0.1 s – 999 h.

1. Change the mode to DISP•Time

2. Enter the required dispensing time with the ▲▼ keys

➡ The display accelerates when pressure on the ▲ or ▼ key is maintained.

Seconds: 0.1"– 899.9" (in 0.1 s steps)

Minutes: 15' – 899' (in 1 min steps)

Hours: 15 h – 999 h (in 1 h steps)

3. Start pump with RUN/STOP

➡ With the ▲▼ keys the dispensing time can be changed while the dispensing cycle is in process.

➡ During the dispensing process the flow rate (PUMP • Flow rate) and the total volume can be re-trieved by pressing the MODE-key. The flow rate can be changed during the dispensing process by using the ▲▼ keys.

➡ Very short dispensing times can result in non-repeatable dispensing volumes.

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
Please note

For entering the number of dispensing cycles  
see Page 13.

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## Please note

For entering the number of dispensing cycles see Page 13.

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## Intermittent dispensing (by time)

### Intermittent dispensing by time with a pre-set pause time

1. Go to mode PAUSE • Time  
Enter the pause time (between 0.1 s – 999 h) with 56 keys.
2. Press MODE key until both LEDs for PAUSE and DISP•Time light up simultaneously.  
➡ Enter the dispensing time (0.1 s – 999 h) with ▲▼ keys.
3. Start with the RUN/STOP key  
➡ The pump dispenses now at intervals.  
➡ Active mode blinks.
4. Pushing the RUN/STOP key again interrupts the dispensing process or pause period.  
➡ The display blinks. The user can now proceed as follows:  
➡ By pressing the RUN/STOP key again, the dispensing or pause period is continued.  
➡ By pressing the RESET-key the pump is stopped for good.  
➡ With the ▲▼ keys the dispensing or pause time can be changed even during the dispensing process. The newly entered set point appears shortly on the display.  
➡ During the dispensing process the flow rate can be retrieved with the MODE key and adjusted with the ▲▼ keys.

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### Please note

With RUN/STOP a dispensing cycle can be interrupted and resumed.  
With RESET an interrupted dispensing cycle can be definitely stopped.

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## Number of dispensing cycles

The number of dispensing cycles can be entered when dispensing at intervals (by time and volume).  
See basic settings on Page 13.

- ➡ During the pauses the display shows the remaining number of dispensing cycles.
- ➡ During long pauses, the remaining interval time can be checked with the CAL/MAX key and adjusted by using the ▲▼ keys.
- ➡ At any time, the dispensing process can be interrupted by pushing the RUN/STOP-key. When resuming the process, the pump continues with the remaining number of dispensing cycles.

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Drip-free dispensing

With programmable roller back-steps (1 – 100 steps), see basic settings bStP (backsteps) Page 13.

Roller back-steps are active in the following modes:

- a. Dispensing by volume (see Page 20)
- b. Dispensing a volume within a pre-set time (see Page 24)
- c. Intermittent dispensing by volume (see Page 25)

⚠ Reject the first dispensing step. Afterwards, the roller back-steps are compensated.

Deleting roller piston strokes  
See basic settings bStP, Page 13.

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
## Please note

The analog interface is integrated in the port 'RS232 OUT' and can be used for remote control (e.g. by footswitch).

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## Analog interface

### Pin 4: +5VDC

For connecting an external control unit (current <100mA).

### Pin 5: GND

Point of reference for all signals described.

### Pin 6: start

Connection for an external trigger device (footswitch, etc.) For operating modes, refer to Basic settings "Footswitch" Page 13.

### Pin 8: busy

Informs about pump status:

Pumpe in operation = output +5V  
Pumpe in stand-by = output on ground (The output resistance is 4.7 kΩ).

### Pin 9: speed OUT

A square-wave signal with a frequency proportional to the speed (rpm) is at your disposal:

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
## Please note

When switching on the REGLO *Digital* pump, the serial interface will only respond after a delay of 5 seconds.

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## Serial interface

### RS232 IN (female)

A 9-pin D-socket is available on the rear panel of the pump.

### Setting the pump address

Each command string begins with the pump address (1 – 8). It is factory-set at 1 and can be changed with @. The address will be stored permanently (even after the pump has been switched off). The assignment of the address enables the user to control up to 8 ISMATEC® drives via one interface.

### Structure of the commands

The address is followed by a character. Some commands have an additional parameter which always consists of 4 or 5 figures.

The command string **is completed by** a carriage return (ASCII 13). The pump confirms most of the commands with an asterix \*. Yes/No inquiries are answered by + (yes) or – (no). Multi-digit replies are concluded by (cr) ASCII 13 and (lf) ASCII 10.

**Incorrect** command strings are answered by #.

**Numerical values** are confirmed as 3 to 5-digit figures. Four of the five digits are numerals, one digit is either a decimal point or a preceding blank space.

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
### RS232 OUT (male)

#### Pin 2: RS232 Rx

Is used for connecting additional pumps together with Pin 3/5.

#### Pin 3: RS232 Tx

Is used for connecting additional pumps together with Pin 2/5.

#### Pin 5: GND

Point of reference for all signals described.

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Each pump requires an individual address  
(see Page 31).

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## Cascading several pumps

Providing that an appropriate software is available, up to 8 ISMATEC® pumps can be controlled via one PC.

In order to establish a cascade you need to connect the pumps with the connecting cable (AG0013) as shown on the left.

The same cable AG0013 serves to connect the first pump to the PC, provided that it has a 9-pin RS232 connector.

If the PC is equipped with a 25-pin connector, you need an additional adapter XC0009.

➡ **Each pump must be allocated its own individual address.**

For this procedure, you first have to connect each pump in turn to the RS232 interface of the PC. Two or more connected pumps must not have the same address. Only when all the pumps have been assigned their addresses can you cascade them as shown on the left.

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Please note

The values indicated in this table are only approximate and determined as follows: ml/min, per channel, with water and Tygon® tubing, without differential pressure.

Order No. = Tygon® ST (R3603)

We will be pleased to help the user to select the tubing material. However, the user has the final responsibility for the selection of the correct tubing material. For reproducible results we recommend you use tubing from ISMATEC®.

Please ask for the detailed tubing documentation.

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
Bestell-Nr. AG0013

## Programming Software

### LabVIEW® from National Instruments®

LabVIEW® is a Windows® software for creating flexible and scalable test, measurement, and control applications in connection with various laboratory instruments.  
<http://www.ni.com>

### LabVIEW® - Driver

The current drivers for your ISMATEC® pump can be downloaded from our web-site: <http://www.ismatec.com>

### Interface card Moxa C168

including a 25 pin multiple cable (male) for controlling up to 8 laboratory instruments

Order No. IM0015C168

### Adaptor 25 F / 9M

25 pins, female / 9 pins, male for IM0015C168 (1 adapter for each controlled instrument).

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### Footswitch

Order No. ISM894

This footswitch serves as a Start/Stop device. It is very useful when using the pump as a dispenser for filling tubes, bottles, etc. Both hands are free for handling bottles etc.

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

MS/CA Anpresshebel  
MS/CA pressure lever  
MS/CA levier de pression

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## Spare cassettes

Cassette MS/CA Click'n'go  
Material **POM-C**  
Order No. **IS3510A**

Cassette MS/CA pressure lever  
Material **POM-C**  
Order No. **IS0649A**

Cassette MS/CA pressure lever  
Material **PVDF**  
Order No. **IS3629A**

➡ The cassettes with pressure lever are available on request. This type of cassette may provides better results both for the Tygon MH tubing and at elevated differential pressure conditions.

➡ The cassettes allow the user to insert tubes with different diameters and materials on the same roller-head.

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### Materials POM-C

Polyoxymethylene-Copolymer

- ▶ good chemical resistance to many organic solvents and strong alkaline chemicals
- ▶ is affected by strong acids and oxidizing substances
- ▶ UV-stabilized and stable up to temperatures of 80°C/176°F (dry, continuous use) or 136°C/277°F (dry, for a short time)

### Materials PVDF

Polyvinylidene fluoride

- ▶ very good chemical resistance to acids and most aliphatic, aromatic and chlorinated solvents
- ▶ not suitable for long contact with esters, ketones, amines and strong alkaline chemicals
- ▶ stable both to UV radiation and temperatures up to 110°C/230°F (continuous use) or 142°C/288°F (for a short time)

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## Disposal

Please retain packaging materials until the product warranty ends. Afterwards please discard packaging materials in an environmentally-friendly manner according to local regulations.

Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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## Service and repairs

Provided the REGLO Digital is operated properly and in compliance with this manual, the tubing is the only part that is subject to wear and tear.

For repairs send the complete pump with detailed description of the failure to your ISMATEC® distributor. Please use the original ISMATEC® packing or a packet of equal quality.

## Replacement parts

For repairs you intend to carry out by yourself (after the warranty period), you can ask from your ISMATEC® distributor:

- Parts lists
- Replacement parts for pump drive
- Wiring diagrams

Please give information on defect, date of purchase, serial No., and model.

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## Please note

We also recommend you observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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Learn more  
about the  
Ismatec® product  
line by visiting:  
[www.ismatec.com](http://www.ismatec.com)

**For ordering and technical support, please contact:**

**North America** [sales@ismatec.com](mailto:sales@ismatec.com) | 1-800-323-4340 | 1-847-549-7600

**Europe** [sales.europe@ismatec.com](mailto:sales.europe@ismatec.com) | +49 (0) 9377 9203-0